

adapt-a-rack™

Flexible multi-tube rack

Each slot adapts to accommodate 12 mm–30 mm (5 mL–50 mL) tubes sizes and holds them in a firm upright position

- Made from POM
- Racks can be linked together for additional holding capacity
- Open sides allow for label viewing
- Flat base to hold tubes upright
- Linked racks should be supported underneath when carrying
- Available in solid or a combination of colors for easy coding identification

- **Autoclavable**

Part No.	Tubes	Wells	Wells Dia.	L x W x H	Pkg.	Rack #1 Color	Rack #2 Color
HS120184	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Blue/Green	Green/Blue
HS120185	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Blue/Yellow	Yellow/Blue
HS120186	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Pink/White	White/Pink
HS120187	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Blue/Blue	Blue/Blue
HS120188	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	White/White	White/White
HS120189	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Pink/Pink	Pink/Pink
HS120190	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Yellow/Yellow	Yellow/Yellow
HS120191	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Green/Green	Green/Green



Available in solid or a combination of colors



Holds **different size tubes** in each slot



Universal Tube Rack


Ideal for holding a wide range of tubes

- Made from polypropylene
- Holds up to 12 similar size tubes
- Imprinted references for easy identification
- Side handles for easy transport
- **Not recommended for autoclaving**



Adjustable slot adapts to different diameter tubes

Part No.	Tubes	Wells	Well Dia.	L x W x H	Pkg.	Color
HS24330A	15-50 mL	12	16-30 mm	7.6 x 5.5 x 2 in (19.4 x 14 x 5.1 cm)	1/ea	Blue/Purple ■ ■



Hold and pour feature



grab rack™

Rack features silicon grips to hold tubes in place

- ABS rack with silicon gripping fingers
- Turn the rack over and the tubes will not fall out
- No assembly required
- **Not recommended for autoclaving**

Part No.	Tubes	Wells	Well Dia.	Dimensions L x W x H	Pkg.	Color
HS24371A	5-10 mL	50	16 mm	10 x 5 x 2 in (25.9 x 12.7 x 4.9 cm)	1/ea	White <input type="checkbox"/>

Conical Centrifuge


Available in natural and pure black for light sensitive samples

- Attached flat caps for easy one handed use
- 1.5 ml versions have a frosted writing area and a syringe pierceable cap
- 100% virgin polypropylene
- **Disposable**

Certified DNase, RNase and Pyrogen free

Large capacity tubes – 15 and 50 mL

- Lot-to-lot traceability with easy tear packaging
- Writing area and flat polyethylene lids for easy labelling
- Imprinted graduations and a self-standing version (only in clear bulk pack)
- Reusable card used for the racked versions, making them more environmentally friendly



	Part No.	Vol.	Dimensions	Graduations	Max RCF Rating	Pkg.		Color
a	HS4426B	15 mL	OD 17 mm Length 120 mm	0.5 mL	11,000 x g	Bulk wrap sterile	50/bag, 10 bags/case	Natural <input type="checkbox"/>
	HS4426R	15 mL	OD 17 mm Length 120 mm	0.5 mL	11,000 x g	Racked sterile	50/rack, 10 racks/case	Natural <input type="checkbox"/>
b	HS4428	15 mL	OD 17 mm Length 120 mm	0.5 mL	9,400 x g	Bulk wrap sterile	50/bag, 10 bags/case	Black <input checked="" type="checkbox"/>
	HS4427B	50 mL	OD 29 mm Length 114 mm	5.0 mL	9,500 x g	Bulk wrap sterile	25/bag, 20 bags/case	Natural <input type="checkbox"/>
c	HS4427R	50 mL	OD 29 mm Length 114 mm	5.0 mL	9,500 x g	Racked sterile	25/rack, 20 racks/case	Natural <input type="checkbox"/>
	HS4429	50 mL	OD 29 mm + 0.3 mm Length 114 mm	5.0 mL	9,400 x g	Bulk wrap sterile	25/bag, 20 bags/case	Black <input checked="" type="checkbox"/>
e	HS4427S	50 mL	OD 29 mm + 0.3 mm Length 114 mm	5.0 mL	9,500 x g	Bulk wrap sterile skirted	25/bag, 20 bags/case	Natural <input type="checkbox"/>