

# Blood Collection Tray

Light weight, heavy on organization

## Blood Collection Tray

Organize your phlebotomy, I.V. start or catheter insertion supplies with this spacious, robust and easy to clean carry tray.

- Lightweight easy-to-clean abs construction
- Customize the compartments to suit your personal needs with the dividers provided.
- Comfort grip handle tested to hold up to 10 kg
- Comes complete with a 17 mm well 60 place tube rack

| Part No.        | L x W x H                                      | Handle Height    | Tube Rack  |
|-----------------|--|------------------|------------|
| <b>HS120261</b> | 16.4 x 12.2 x 3.3 in<br>(41.9 x 31.1 x 8.5 cm) | 9.4 in (24.1 cm) | 17 mm rack |



Comes complete with a 17 mm well, 60 place, tube rack





## Sharps Containers

### Easy to use and transport

- Provides safe disposal of syringes, needles and other small sharps
- Stackable to reduce storage space before use
- 1.0 L capacity (quart) with a clear lid to view fill level
- Slide access lid for temporary and final closures for safety and security
- Built-in needle removal port, prevents direct contact with sharps
- Compatible with the droplet® blood collection tray, totemax™ and blood collection tray HS120261

| Part No.        | L x W x H                                  | Pkg.  | Color  |
|-----------------|--|-------|--|
| <b>HS120177</b> | 4.5 x 4.5 x 7.5 in (11.5 x 11.5 x 19.2 cm) | 18/pk | Red <span style="color: red;">■</span>       |
| <b>HS120178</b> | 4.5 x 4.5 x 7.5 in (11.5 x 11.5 x 19.2 cm) | 18/pk | Yellow <span style="color: yellow;">■</span> |



## x-tourn™ tourniquets

### Hi-viz orange alert for safer practice

High-quality disposable tourniquets eliminate roping or curling when used.

Latex free, no rubber proteins, resins or fatty acids that can produce hypersensitive reactions.

- Made from and smooth polyisoprene
- Less abrasive than other brands, uniform thickness prevents shredding and breaking

### - Disposable

| Part No.        | L x W                       | Pkg.    | Color  |
|-----------------|-----------------------------|---------|--|
| <b>HS120108</b> | 18 x 1 in (45.72 x 2.54 cm) | 100/pk  | Green <span style="color: green;">■</span>   |
| <b>HS120109</b> | 18 x 1 in (45.72 x 2.54 cm) | 25/roll | Green <span style="color: green;">■</span>   |
| <b>HS120110</b> | 18 x 1 in (45.72 x 2.54 cm) | 250/pk  | Orange <span style="color: orange;">■</span> |
| <b>HS120111</b> | 18 x 1 in (45.72 x 2.54 cm) | 25/roll | Orange <span style="color: orange;">■</span> |



Available in bulk or on a perforated roll