

Microtube Boxes & Racks

Safe microtube usage and storage



Microtube tube racks

- Made from reinforced nylon
- Superior chemical and heat resistance
- Handles provide easy and safe transportation
- Stackable two formats available for smaller and larger requirements
- Temperature range -40°C to +140°C

- **Autoclavable**

| Part No. | Tubes | Wells | Well Dia. | L x W x H | Pkg. | Color |
|-----------------|------------|-------|-----------|--------------------------------------|------|--|
| HS120241 | 1.5/2.0 mL | 24 | 11.1 mm | 5.6 x 3.5 2 in (14.3 x 9 x 5.2 cm) | 8/pk | Blue ■ |
| HS120242 | 1.5/2.0 mL | 24 | 11.1 mm | 5.6 x 3.5 2 in (14.3 x 9 x 5.2 cm) | 8/pk | White <input type="checkbox"/> |
| HS120243 | 1.5/2.0 mL | 96 | 11.1 mm | 8.2 x 5.1 x 2.1 (21 x 13.1 x 5.4 cm) | 4/pk | Blue ■ |
| HS120244 | 1.5/2.0 mL | 96 | 11.1 mm | 8.2 x 5.1 x 2.1 (21 x 13.1 x 5.4 cm) | 4/pk | White <input type="checkbox"/> |

Microtube racks 80-well

Popular fraction-collector-plate style rack

- Made of polypropylene
- Easily fill and transfer 80 tubes in compact arrangement
- Compact arrangement (5 x 16 array)
- Molded-in coordinates to aid organization

- **Autoclavable**

| Part No. | Tubes | Wells | Well Dia. | L x W x H | Pkg. | Color |
|-----------------|------------|-------|-----------|--|------|--|
| HS29025A | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | Natural <input type="checkbox"/> |
| HS29025B | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | Blue ■ |
| HS29025C | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | Green ■ |
| HS29025D | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | Pink ■ |
| HS29025E | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | Yellow ■ |
| HS29025F | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | Orange ■ |
| HS29025G | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | ■ ■ ■ ■ ■ * |
| HS29025J | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | Red ■ |

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



Cardboard micro tube boxes

Safe micro tube storage

Low profile boxes are available in two sizes to accommodate standard 0.5 mL and 0.2 mL micro tubes and vessels of a similar size.

- Dividers included with the boxes
- Store in mechanical freezers or vapor phase LN2*
- Will fit into most standard freezer racking

| Part No. | Tubes | Wells | Well Dia. | L x W x H | Pkg. | Color |
|-----------------|--------|-------|-----------|--|-------|--------------------------------|
| HS120220 | 0.2 mL | 169 | 9.5 mm | 5.2 x 5.2 x 1.1 in (13.3 x 13.3 x 2.9 cm) | 10/pk | White <input type="checkbox"/> |
| HS120221 | 0.5 mL | 100 | 12 mm | 5.2 x 5.2 x 1.6 in (13.3 x 13.3 x 4.1 cm) | 10/pk | White <input type="checkbox"/> |

*Must use tubes capable of being stored at low temperatures



Cardboard cryogenic vial boxes and partitions

For storage of small samples and tubes

- Store in LN2, or in mechanical freezers
- Colors allow for easy identification and sorting purposes
- Bottom vent holes to allow quicker circulation and removal of cold air, vapor or LN2
- All box footprints are 13.3 x 13.3 cm, select from two different heights

| Part No. | Description | Tubes | Wells | Wells dia. | L x W x H | Pkg. | Color |
|-----------------|---------------------|------------|-------|------------|---|-------|---------------------------------|
| HS2860A | Vial box and lid | 1.5/2.0 mL | — | — | 5.2 x 5.2 x 2 in (13.3 x 13.3 x 5 cm) | 12/pk | White <input type="checkbox"/> |
| HS2860B | Vial box and lid | 1.5/2.0 mL | — | — | 5.2 x 5.2 x 3 in (13.3 x 13.3 x 7.5 cm) | 12/pk | White <input type="checkbox"/> |
| HS2860C | Vial box and lid | 1.5/2.0 mL | — | — | 5.2 x 5.2 x 2 in (13.3 x 13.3 x 5 cm) | 12/pk | Red, Blue, Yellow* |
| HS2860CB | Vial box and lid | 1.5/2.0 mL | — | — | 5.2 x 5.2 x 2 in (13.3 x 13.3 x 5 cm) | 12/pk | Blue <input type="checkbox"/> |
| HS2860CR | Vial box and lid | 1.5/2.0 mL | — | — | 5.2 x 5.2 x 2 in (13.3 x 13.3 x 5 cm) | 12/pk | Red <input type="checkbox"/> |
| HS2860CY | Vial box and lid | 1.5/2.0 mL | — | — | 5.2 x 5.2 x 2 in (13.3 x 13.3 x 5 cm) | 12/pk | Yellow <input type="checkbox"/> |
| HS2860E | Partition (10 x 10) | 1.5/2.0 mL | 100 | 11.5 mm | — | 12/pk | |
| HS2860H | Partition (9 x 9) | 1.5/2.0 mL | 81 | 12 mm | — | 12/pk | |
| HS2860F | Partition (8 x 8) | 1.5/2.0 mL | 64 | 15 mm | — | 12/pk | |
| HS2860G | Partition (7 x 7) | 1.5/2.0 mL | 49 | 17.5 mm | — | 12/pk | |

*Assorted colors (Red, Blue, Yellow)

Partitions must be ordered separately