


Fully assembled, polypropylene mega racks are great for high throughput labs requiring large sample storage


- Sturdy, durable material for easy and secure transport of clinical and research tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Reusable, easy to clean and Autoclavable
- Place different test groupings into one rack eliminating the need for multiple racks
- Stackable for staging and saving space



120114

Mega Racks® 10-13 mm - Holds 5-7 mL sample tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24322A	Blue ■	216	13	13.3 x 9.3 x 3.6	33.7 x 23.5 x 9.1
HS24322E	White □				
120112	White □	432	13	26.5 x 9.3 x 3.5	67.3 x 23.5 x 9.0
120113	Blue ■				



HS24323A & HS24323E


Mega Racks® 13-16 mm - Holds 10 mL sample tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24323A	Blue ■	120	16	13.3 x 10.5 x 3.6	33.7 x 26.7 x 9.1
HS24323E	White □				
120114	White □	240	16	26.5 x 10.5 x 3.5	67.3 x 26.7 x 9.0
120115	Blue ■				

Mega Racks® 18-20 mm - Holds 15 mL conical/swab tubes and 20 mm round bottom tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120834	White □	120	20	14.9 x 11.9 x 3.6	37.7 x 30.1 x 9.1
120835	Blue ■				

Ideal footprint
for use with programmable tube filling machines



120834 & 120835

