


Three Scales: Honey Moisture, Brix, Baume

Simple to operate and provides quick, accurate measurements for honey sampling

Key Features

- Large, Sharp Display Scales with high contrast
- High-quality optical glass material for consistent, accurate readings
- Durable, industrial-grade design with non-roll stand feature to keep unit secure on flat surfaces
- Tactical, non-slip Rubber Grip on body to block transmission of hand heat, reducing temperature effects
- Prism encased in all-metal housing, allowing sample and prism to reach temperature equilibrium quickly
- Five-Year Warranty



Item No.	Model	Brix range	Water Range	Baume Range	ATC
43006	HMX-1	58.0-92.0%	13-27%	38-43°	Yes

Replacement parts	
43RP-008	Soft Fabric Carrying Case for all Handheld Refractometers
45000-OS025	Refractometer Standard, Oil-based, 60.00% Brix / 1.4420 nD (20°C), 7.4 mL/0.25 oz
43000-DLHMX1	Replacement Daylight Plate Kit, for Model HMX-1 refractometer

Specifications	
Range	Brix: 58.0-92.0% / Water: 13-27% / Baume: 38-43°
Resolution	Brix: 0.5% / Water: approximately 1% / Baume: 0.5°
Accuracy	Brix: ±0.5% / Water: ±1% / Baume: ±0.5°
ATC Range	10-30°C (50-86°F)
Dimensions	1.6 x 1.6 x 5.8 in (4 x 4 x 14.8 cm)
Weight	220 g (7.8 oz)
Max Sample Temp.	90° C (194° F)

Package Includes	
Honey Refractometer	
Vinyl storage/carrying case	
One transfer pipet	
Calibration screwdriver	
Operation manual	