

DuraPorter® XL High-Capacity Sample Transport Tote

Your reliable solution for the safe and convenient transportation of laboratory samples, instruments, and valuable items. Crafted to meet the stringent demands of clinical and research laboratories, this oversized transport container ensures the protection and preservation of your cargo.

Built to Last and Designed with the User in Mind

- Sturdy and spacious, ideal for transporting a large volume of samples
- Transparent design helps comply with OSHA 29 CFR 1910.1030 standards without having to re-label containers and enabling the user to visibly check sample integrity and meets UN3373 secondary packaging requirements
- Equipped with articulated handles featuring silicone grips for comfortable and secure handling



Emphasizing Security and Safety

- A silicone seal ensures a tight enclosure, safeguarding samples from most environmental elements
- Dual latches with security tag feature guarantee a secure closure, providing tamper-evidence and ensuring the integrity and authenticity of your samples



Exceptional Versatility

- Includes a MultiTube Rack Insert to accommodate various tube sizes. *Additionally, two other insert options for blood tubes up to 10 mL and general research tubes up to 15 mL are available separately*
- The bottom compartment is generously sized, allowing for carrying most tube racks, instruments or supplies according to your specific requirements

Soft Grip Handle
for comfortable handling

Translucent Design
allows for easy content
identification

Removable Lid
for ease of loading and
unloading



Silicone Seal
for secure content enclosure

4-Point Handle Security
ensures durable transport

Image shown with 120961 Small Tube Rack Insert (not included)

Dual Latches
with tag out feature for
added security

Available in 3 Colors
for easy user or sample
identification and organization



Includes One 120962 MultiTube Rack Insert

Item No.	Color	Dimensions (not including handle) in	Dimensions (not including handle) cm	UOM
120926	Blue	16 x 16.25 x 7	40.64 x 41.28 x 17.78	1 ea
120927	Red	16 x 16.25 x 7	40.64 x 41.28 x 17.78	1 ea
120928	Yellow	16 x 16.25 x 7	40.64 x 41.28 x 17.78	1 ea

Duraporter® XL
 MultiTube Rack Insert



- Designed to fit DuraPorter® XL High-Capacity Sample Transport Tote
- Thermoformed PETG
- Tube Rack Inserts can be used separately
- Lid removes completely for easy cleaning
- Insert deep enough to be filled with ice

Item No.	Tubes	Cut-out Diameter	UOM
120960	15 mL, Large blood tube	18 mm	Pack of 3
120961	Small blood tubes	13 mm	Pack of 3
120962	50 mL, 15 mL, Small & Large blood tubes	30 mm (center), 18 mm (large), 13 mm (small)	Pack of 3



120960

Duraporter® XL Large Tube Insert

Holds 100 tubes with diameter up to 18 mm
 Most 10 mL blood tubes & 15 mL conical tubes
 Array 10 x 10
 Pack of 3



120961

Duraporter® XL Small Tube Insert

Holds 168 tubes with diameter up to 13 mm
 Most 3.5 - 7.5 mL blood tubes
 Array 12 x 14
 Pack of 3



120962

Duraporter® XL MultiTube Rack Insert

Holds tube with various diameters: 13, 18 & 30 mm
 Most 15 mL & 50 mL conical tubes
 and 3.5 mL - 10 mL blood tubes
 Pack of 3