

Lightweight, roomy, and easy-to-clean blood collection tray

Designed with phlebotomists in mind!

Keeping a clean and organized tray has never been easier.

Fully equipped with tray inserts, tube rack, supply drawer, side pockets, and a ratcheting handle with soft rubber grip that rotates out of the way for easy access.

Tray inserts come in two configurations based upon the size of tube rack:

- Choice of 36-Place 13 mm or 16 mm Heathrow OneRack® test tube rack
- Features an extra-deep compartment for a sharps container
- Made of easy to clean ABS material
- Replacement tray inserts available



Droplet™ Configurations

Part No.	Description	L x W x H	Handle Height
HS2200A	Droplet™ blood collection tray with 2 inserts (style A) and a 36-place, 13 mm test tube rack	16.8 x 10.4 x 6.7 in (42.9 x 26.5 x 17.2 cm)	11 in (28 cm)
HS2200B	Droplet™ blood collection tray with 2 inserts (style B) and a 36-place, 16 mm test tube rack	16.8 x 10.4 x 6.7 in (42.9 x 26.5 x 17.2 cm)	11 in (28 cm)

Droplet™ Inserts

Part No.	Description	L x W
HS2201A	Droplet™ blood collection tray insert A for 13 mm tube rack	14.2 x 8.2 in (36 x 21 cm)
HS2201B	Droplet™ blood collection tray insert B for 16 mm tube rack	14.2 x 8.2 in (36 x 21 cm)



US Patent No.D904026