



Sustainable, Safe Solutions

Reusable lab mats are an essential lab solution to help keep benchtops organized and clean

They help facilitate workflow processes and are ideal for every type of laboratory

Save Money & Eliminate Waste

 Since the lab mats are reusable, they reduce cost and waste making it a more environmentally friendly option over the alternative paper mats

Adds Stability & Safety

- These mats offer a large working surface that is anti-skid which keeps delicate glassware and expensive equipment in place
- The mat's surface is also noise dampening, reducing sound from the vibration of common benchtop equipment
- Made from durable food grade silicone, they are chemically resistance and have a raised edge to contain chemical spills

Easy to Clean & Store

- When the mats are done being used, they can be easily cleaned using any standard disinfectant by simply wiping it, or putting it in a dishwasher and it's ready for repeated use



Workstation Lab Mat Safely facilitate workflow processes

- Protect counter tops from hot items up to 200°C (392°F)
- Lip around the edge can hold approx. 250 mL of liquid
- Small footprint, 12 x 14.9 x 0.1 in (30.5 x 37.9 x 0.2 cm)





Gridded Surface

for tube staging, preventing movement of objects



Side One: combination gridded and smooth surface
Side Two: smooth surface only





Item No.	Color	L x W x H in	L x W x H cm
120747	Blue	12 x 14.9 x 0.1	30.5 x 37.9 x 0.2
120748	Purple	12 X 14.9 X U.1	30.3 X 37.9 X 0.2

Silicone Lab Mat

Large working space

The lab mat can be used with items such as centrifuges, vortex mixers, hot plates, stirrers, pipettes and so much more $\frac{1}{2}$

- Mat features designated magnetic area to keep stir bars in place
- Protect countertops from hot items up to 200°C (392°F)



Raised Edges

enable liquid retention of approximately 300 mL of fluid



Item No.	Side #1 Color	Side #2 Color	L x W x H in	L x W x H cm
120506	Yellow	Grey/Blue		
120507	Purple	Purple/Grey	14 x 23.5 x 0.1	35 x 60 x 0.2
120500	Rlue	White/Blue		




Reversible, with two color/design options available for color coding or to add color to work spaces

Large Working Surface

Solid Color on both sides!

Design includes internal metal inserts for retention of magnetic objects, such as stir bars to help keep items organized





Lip design around all edges enables liquid retention of approximately 300mL of fluid





