

Microtube Racks

Durable, Simple to Use

- Made from reinforced nylon offering superior chemical and heat-resistance over most other plastic based racks
- Handles provide easy and safe transportation
- Stable design can be used in water baths
- Alphanumeric identification for quick location and retrieval of tubes



| Item No. | Color | Tubes mL | Wells | Well Dia. mm | Packaging | L x W x H in | L x W x H cm |
|----------|-------|-------------|-------|-----------------|-----------|-----------------|-----------------|
| 120241 | Blue | ■ 1.5/2.0 | 24 | 11.1 | 8 | 5.6 x 3.5 x 2 | 14.3 x 9 x 5.2 |
| 120243 | | | 96 | | 4 | 8.3 x 5.2 x 2.1 | 21 x 13.1 x 5.4 |

Microtube Racks 80-well

Popular fraction-collector-plate style rack

Compact 5x16 array makes filling and transfer easy

- Frosted labeling area on bottom of rack
- Embossed index numbers for quick sample ID and retrieval
- Made of polypropylene



Pack of 5

| Item No. | Color | Part No. | Color |
|----------|-----------|----------|---------------------|
| HS29025A | Natural □ | HS29025E | Yellow ■ |
| HS29025B | Blue ■ | HS29025F | Orange ■ |
| HS29025C | Green ■ | HS29025G | Assorted* ■ ■ ■ ■ ■ |
| HS29025D | Pink ■ | HS29025J | Red ■ |

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|------------------|
| 1.5/2.0 | 80 | 11 | 8.9 x 2.6 x 1 | 22.5 x 6.7 x 2.8 |





Microtube storage boxes

Safe microtube storage

Low profile boxes are available in two sizes to accommodate standard 0.5 mL and 0.2 mL microtubes and vessels of a similar size

- Partitions included in the boxes
- Will fit into most standard freezer racking



Pack of 10

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|----------|-------|-------------|-------|-----------------|-----------------|-------------------|
| 120220 | White | □ 0.2 | 169 | 9.5 | 5.2 x 5.2 x 1.1 | 13.3 x 13.3 x 2.9 |
| 120221 | White | □ 0.5 | 100 | 12 | 5.2 x 5.2 x 1.6 | 13.3 x 13.3 x 4.1 |

Cryogenic Vial Boxes & Partitions

For storage of small samples and tubes

- Made from cardboard
- Colors allow for easy identification and sorting purposes
- All box footprints are 13.3 x 13.3 cm, select from two different heights



Vent Holes

allow faster cooling and removal of LN₂

Cardboard Cryogenic Vial Boxes

Pack of 12

| Item No. | Color | Tubes mL | L x W x H in | L x W x H cm |
|----------|----------|-------------|-----------------|-------------------|
| HS2860A | White | □ | 5.2 x 5.2 x 2 | 13.3 x 13.3 x 5 |
| HS2860B | White | □ | 5.2 x 5.2 x 3 | 13.3 x 13.3 x 7.5 |
| HS2860C | Assorted | ■ ■ ■ | | |
| HS2860CB | Blue | ■ | | |
| HS2860CR | Red | ■ | 5.2 x 5.2 x 2 | 13.3 x 13.3 x 5 |
| HS2860CY | Yellow | ■ | | |

Partitions

Partitions must be ordered separately

Pack of 12

| Item No. | Array | Tubes mL | Wells | Well Dia. mm |
|----------|---------|-------------|-------|-----------------|
| HS2860E | 10 x 10 | 1.5/2.0 | 100 | 11.5 |
| HS2860H | 9 x 9 | 1.5/2.0 | 81 | 12 |
| HS2860F | 8 x 8 | 1.5/2.0 | 64 | 15 |
| HS2860G | 7 x 7 | 1.5/2.0 | 49 | 17.5 |