



# Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (76 x 26 mm and 3 x 1 in)

## Organize & Identify

- Lining and slots are numbered to correspond to the slide inventory sheet
- Organize by Project and Date with pre-molded labeling area
- Available in multiple colors for project coding

## Durable & Safe

- ABS plastic boxes available with either cork or foam lining
- Designed to securely store slides in grooved slots
- Molded feet secure stacking for easy transport and storage



## 25-Place Slide Box

Easily access slides with friction- t cover



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15989A	Blue	Cork	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5
HS15989B	Green				
HS15989C	Red				
HS15989E	White	Foam	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5
HS15989M	Blue				
HS15989P	Red				
HS15989R	White				

## 50-Place Slide Box

Thumb-latch lock mechanism



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15996A	Blue	Cork	Thumb-Latch	8.2 x 3.4 x 1.4	20.8 x 8.6 x 3.5
HS15996B	Red				
HS15996C	White				

## 100-Place Slide Box

Closes securely with a rust-resistant nickel-plated clasp and pin hinges

Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15994A	Blue	Cork	Nickel-plated Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4
HS15994B	Green				
HS15994C	Red				
HS15994D	Yellow				
HS15994E	White				
HS15994F	Grey				
HS15994G	Black	Foam	Nickel-plated Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4
HS15994M	Blue				
HS15994P	Red				
HS15994R	White				

## 100-Place Tall Slide Box

Additional height for 3 x 2 in slides  
Stackable with 100-Place Premium Slide Boxes

Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15979A	Blue	Cork	Nickel-plated Clasp	8.2 x 6.9 x 2.4	20.8 x 17.5 x 6