

Packed with features for comfortable and accurate pipetting

The powerful pump offers rapid, yet sensitive aspiration and dispensing, has adjustable valves and is microprocessor controlled for constant reliability and accuracy

- Long-life li-ion battery, with quick charge—up to 20 hours use
- Toggle operation by pressing the large mode button
- An adapter can be fitted to the nose cone to stabilize small volume and thinner pipettes
- New USB connector ensures the appropriate charger is used with the unit
- Low-battery indicator
- Can be recharged when in use



Rechargeable li-ion battery (replaceable)



Nose cone assembly only

Item No.	Volume Range mL	L x W x H in	L x W x H cm
HS3000	1–100	9 x 8.5 x 5	22.9 x 21.6 x 12.7

Replacement parts

HS5000RFFIL	RF3000® Replacement Filters, Pack of 5
120669	RF3000® Universal Power Adaptor 5VDC Micro USB Plug
HS300502	RF3000® Li-Ion Rechargeable Battery, One year warranty on battery included
120487	RF3000® Nose Cones Assorted* ■ ■ ■ ■ ■ Set of 5
HS5000RFRFA	RF3000® Replacement Filter Adapter
HS5000RFRP	RF3000® Replacement Pipette Silicone Adapter

*Assorted colors (Black, Blue, Grey, Red and Yellow)

Specifications

Power	Micro-USB connector Rechargeable li-ion battery (replaceable)
Filter	PTFE/25mm dia/0.45 µm filter
Pipettes	1–100 mL, glass and plastic
Pipette Nozzle	Polypropylene

Package Includes

RF3000® Pipette Controller
Charging Stand and Wall-Mountable/Magnetic Holder
Nozzle Adapter for use with Small Pipettes
Li-Ion Rechargeable Battery
Five Nose Cones: Black, Yellow, Red, Grey and Blue
5VDC Micro USB Plug
6 ft (183 cm) Cable

1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



- Soft-Coated Triggers with direction indicators
- Universal Grip for right- or left-handed users
- Ergonomically Balanced in the hand with a loaded pipette
- UV resistant housing

Easy Identification

5 different colored nose cones help differentiate department, task or user; as well as a labeling area on the side of the unit

Adapter for small pipettes

Storage Stand doubles as a charging station



Magnetic & Screw Wall-Mount



3 Modes
Bright LEDs indicate 3 operation modes: high, low, & gravity