



Weighted Stand
for bench top use

Lock Counter Off
to allow use of the marker
without registering more counts

Hand Tally Adapter
allows you to use the eCount
as a digital tally counter

eCount Colony Counter

Combines a multi-function electronic counter with a Sharpie® Pen for marking to prevent missing or double-counting colonies

- Light, ergonomic body is balanced for easy manipulation
- Select between an audible beep or an LED flash to acknowledge count
- Store up to 32 separate counts and display the total of all counts
- Count up/down verification
- Quick and easy to change between different colored markers to check different colony counts in the same view



Battery and Pen not included in warranty



US Patent No. D564,385
EP Registered Design No. 705751-0001

Item No.	L x W in	L x W cm
120000	7.4 x 1.2*	17 x 3*
120228	CR1225 Battery	

*Measurements shown are with Sharpie® Pen inserted

Package Includes

eCount Colony Counter	Hand Tally Adapter
Black Sharpie® Pen	Weighted Stand
Long-life Button Battery	

