

Ideal for microbalance weighing and/or distribution of very small amounts of sampling material

Easy to Handle


- Large tabs make boats easy to use, position and hold, especially with forceps

Superb Sample Recovery

- Inert material enables the boat and contents to be added to solutions to aid 100% sample recovery, where protocols allow
- Anti-static properties prevent sample retention from static electricity
- Small surface area reduces liquid sample evaporation

Micro Weigh Boats

Item No.	Color	Dimensions <i>in</i>	Dimensions <i>mm</i>	UOM
120836	White	.98 x 0.4 x 0.3	25 x 8.9 x 7.5	250/pk



Actual size



Easy to Handle Tabs

Made of Polystyrene
Food Grade Compliant

Holds approx 500 µl



Ideal for very small volume
pipette calibration and
checks and micro sampling

