



The dual purpose scoops by Heathrow Scientific are a unique economical solution that combines the function of a scoop weigh boat, and funnel into one product.

The weighing process has been streamlined, now just sample or retrieve the material in the scoop and place it directly on the scale. This eliminates the need for a weigh boat or weighing paper and saves valuable space in the lab.

- Made from anti-static polystyrene that is waterproof and reduces the chance of materials sticking to the scoop
- Ships and stores flat and when ready to use, is easy to assemble by inserting the locking tab into the hole
- Locking tab feature enables a scoop shape to help retain materials
- Coated to stand up to powders, gels, and semi-solids
- Broad tip design keeps the scoop stable while on the scale and enables materials to be easily scooped and dispensed
- Disposable



Item No.	Color	Size	L x W in	L x W cm	UOM
120656	White	<input type="checkbox"/> Large	5 x 3.5	12.7 x 9	100 ea/pk, 50 pks/cs
120657	White	<input type="checkbox"/> Small	5 x 2.5	12.7 x 6.3	100 ea/pk, 50 pks/cs



Easily assemble by pulling the arrow tab through the hole.

