

**Each well has fingers that flex to accommodate a variety of tube sizes while holding them in a firm upright position**

- Made of POM
- Racks can be linked together for additional capacity
- Open sides allow for label viewing
- Flat base to hold tubes upright
- Available in solid or a combination of colors for easy coding identification



Pack of 2

Item No.	Rack #1 Color	Rack #2 Color
120184	Blue/Green	Green/Blue
120185	Blue/Yellow	Yellow/Blue
120186	Pink/White	White/Pink
120187	Blue/Blue	Blue/Blue
120188	White/White	White/White
120189	Pink/Pink	Pink/Pink
120190	Yellow/Yellow	Yellow/Yellow
120191	Green/Green	Green/Green

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
7-50	4	12-30	7.1 x 2.2 x 3	18.1 x 5.6 x 7.6

Holds different size tubes in each slot

