



Premium Plus Slide Boxes

Modern design for daily use

NEW!





- a.** Slide rail design ensures a secure fit allowing easy placement and removal of slides
- b.** Centrally positioned location numbering grid enables simple slide organization with easy viewing of location numbers
- c.** Entirely made from ABS and Polypropylene makes the boxes robust and eliminates potential rusting of metal parts
- d.** Sliding lock provides secure storage with easily identified lock and unlock symbols. 2 mm zip tie can be used for more secure storage
- e.** Inventory card is removable and replaceable for trouble-free writing on flat surfaces. Inventory cards can be purchased separately in pks/10
- f.** Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement
- g.** EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- h.** Multiple labeling areas, a large area on the lid of the box with “project” and “date” molded in and the back of the index card can be used as the lid is transparent



- i.** Vibrant translucent colors are easy to identify and enable coding options, assorted 5/pk available
- j.** Stylish wave surface texture is easy to hold, increasing safety of sample transport
- k.** Stackable, molded feet allow secure horizontal stacking of multiple units for safe transport or storing

- l.** Flexible storing, the units can be stored vertically or horizontally helping to efficiently use space available
- m.** Unique side storage can accommodate a pen, desiccant or other small items needed for working
- n.** Holds 100 - 26 x 75 mm slides

Patent Pending



NEW!



100-Place Premium Plus Slide Box

| Item No. | Color | Lining | Hinge | L x W x H <i>in</i> | L x W x H <i>cm</i> | UOM |
|----------|-----------|----------|---------------|------------------------|------------------------|---|
| 120576 | Blue | EVA Foam | Sliding Latch | 8.2 x 7.7 x 1.3 | 20.7 x 19.5 x 3.2 | 5/pk 4/pks per case (shrink wrapped packs) |
| 120577 | Orange | | | | | |
| 120578 | Green | | | | | |
| 120579 | Red | | | | | |
| 120592 | Purple | | | | | |
| 120593 | Assorted* | | | | | |

*Assorted colors (Purple, Blue, Green, Red, Orange)

Inventory Card

| Item No. | L x H <i>in</i> | L x H <i>cm</i> | UOM |
|----------|--------------------|--------------------|--|
| 120618 | 8.5 x 11 | 21.6 x 28 | 10/PK (Polybag with Label)50 PK per case |