



Unique pipette holder that gives you the option of using on the bench top or hanging under shelf to help maximize your working space

- Will hold up to 12 single-channel pipettes, or a combo of 3 single- and 3 multichannel pipettes
- Holds most common brands
- Rubber grips prevent scratching pipette
- Weighted base
- Installation hardware included

US Design Patent No. D754,872 & D763,464

Item No.	Color	Bench Top Style W x H in	Bench Top Style W x H cm	Under Shelf Style W x H in	Under Shelf Style W x H cm
120026	Blue				
120056	Lime	13.1 x 7.6	33.2 x 19.3	7.4 x 7.6	18.7 x 19.3
120057	Magenta				

