

Lightweight and easy to use

The RF1000 is a manual pipette controller uniquely designed for one-handed operation. Simply squeeze the main bulb and push the toggle forward to aspirate or back to dispense. Push the blowout bulb to expel any liquid remaining in the pipette.

- Tapered silicone rubber pipette adapter allows a large range of pipette sizes (1-100 mL) Replaceable when worn.
- Accepts common and readily available replacement filters unlike other pipette aids whose filters are custom. This will help reduce total cost of ownership over time.
- Includes one 0.45 µm hydrophobic filter.

Item No.	Description
HS1000RFU	RF1000 Manual Pipette Controller
Replacement parts	
HS5000RFIL	Replacement Filters, pack of 5
120604	Nose Cone with Adapter, Blue
Specifications	
Power:	Manual
Filter:	PTFE/25 mm dia / 0.45 µm filter in PP housing
Pipettes:	1 – 100 mL, glass and plastic
Pipette Cone:	Polypropylene, autoclavable
Dimensions:	180 x 45 x 120 mm



Filter not included
in warranty



Nose cone assembly only

