




## Economical, Reusable & Moisture Resistant Plastic

True North® freezer boxes by Heathrow Scientific, were designed with sustainability in mind and are made from a specially formulated corrugated or thin film polypropylene. True North freezer boxes are moisture resistant, which significantly reduces the chance of mold growth. With superior durability, polypropylene outlasts the cardboard box alternatives that are susceptible to softening and producing dust and can freeze (down to temperatures of  $-80^{\circ}\text{C}$ ) and thaw multiple times without a problem.

When it comes to saving valuable lab space, True North Flatpack Freezer Boxes are supplied “flat-packed”, which eliminates the hassle of stacking and reduces the amount of space needed for storage. When ready to use, assembly is easy and pre-assembled dividers are included with each box. True north freezer boxes come complete in packs of 10, and multiple sizes and colors available.



**Freezer boxes are supplied “flat-packed”**  
Eliminating the hassle of stacking and reduces the amount of space needed for storage





Corrugated Polypropylene 0.2/0.5/1.5-2.0/5 mL

Part No.	Tube	Wells	Well Dia.	L x W x H	Pkg.	Color	
<b>HS120379</b>	0.2 mL	144	7.9 mm	5.2 x 5.1 x 1.2 in (13.3 x 13.0 x 3.0 cm)	10/pk	Natural	<input type="checkbox"/>
<b>HS120380</b>	0.2 mL	144	7.9 mm	5.2 x 5.1 x 1.2 in (13.3 x 13.0 x 3.0 cm)	10/pk	Blue	<input checked="" type="checkbox"/>
<b>HS120381</b>	0.2 mL	144	7.9 mm	5.2 x 5.1 x 1.2 in (13.3 x 13.0 x 3.0 cm)	10/pk	Red	<input checked="" type="checkbox"/>
<b>HS120382</b>	0.2 mL	144	7.9 mm	5.2 x 5.1 x 1.2 in (13.3 x 13.0 x 3.0 cm)	10/pk	Purple	<input checked="" type="checkbox"/>
<b>HS120383</b>	0.5 mL	81	12 mm	5.2 x 5.1 x 1.6 in (13.3 x 13.0 x 4.2 cm)	10/pk	Natural	<input type="checkbox"/>
<b>HS120384</b>	0.5 mL	81	12 mm	5.2 x 5.1 x 1.6 in (13.3 x 13.0 x 4.2 cm)	10/pk	Blue	<input checked="" type="checkbox"/>
<b>HS120385</b>	0.5 mL	81	12 mm	5.2 x 5.1 x 1.6 in (13.3 x 13.0 x 4.2 cm)	10/pk	Red	<input checked="" type="checkbox"/>
<b>HS120386</b>	0.5 mL	81	12 mm	5.2 x 5.1 x 1.6 in (13.3 x 13.0 x 4.2 cm)	10/pk	Purple	<input checked="" type="checkbox"/>
<b>HS120262</b>	1.5/2.0 mL	81	12 mm	5.2 x 5.1 x 2.0 in (13.3 x 13.0 x 5.1 cm)	10/pk	Natural	<input type="checkbox"/>
<b>HS120376</b>	1.5/2.0 mL	81	12 mm	5.2 x 5.1 x 2.0 in (13.3 x 13.0 x 5.1 cm)	10/pk	Blue	<input checked="" type="checkbox"/>
<b>HS120377</b>	1.5/2.0 mL	81	12 mm	5.2 x 5.1 x 2.0 in (13.3 x 13.0 x 5.1 cm)	10/pk	Red	<input checked="" type="checkbox"/>
<b>HS120378</b>	1.5/2.0 mL	81	12 mm	5.2 x 5.1 x 2.0 in (13.3 x 13.0 x 5.1 cm)	10/pk	Purple	<input checked="" type="checkbox"/>
<b>HS120387</b>	5 mL	25	22.6 mm	5.2 x 5.1 x 2.9 in (13.3 x 13.0 x 7.5 cm)	10/pk	Natural	<input type="checkbox"/>
<b>HS120388</b>	5 mL	25	22.6 mm	5.2 x 5.1 x 2.9 in (13.3 x 13.0 x 7.5 cm)	10/pk	Blue	<input checked="" type="checkbox"/>
<b>HS120389</b>	5 mL	25	22.6 mm	5.2 x 5.1 x 2.9 in (13.3 x 13.0 x 7.5 cm)	10/pk	Red	<input checked="" type="checkbox"/>
<b>HS120390</b>	5 mL	25	22.6 mm	5.2 x 5.1 x 2.9 in (13.3 x 13.0 x 7.5 cm)	10/pk	Purple	<input checked="" type="checkbox"/>



Corrugated Polypropylene 15 mL and 50 mL

Part No.	Tube	Wells	Well Dia.	L x W x H	Pkg.	Color	
<b>HS120391</b>	15 mL	36	19.8 mm	5.78 x 5.7 x 4.72 in (14.7 x 14.5 x 12.0 cm)	10/pk	Natural	<input type="checkbox"/>
<b>HS120392</b>	15 mL	36	19.8 mm	5.78 x 5.7 x 4.72 in (14.7 x 14.5 x 12.0 cm)	10/pk	Blue	<input checked="" type="checkbox"/>
<b>HS120393</b>	15 mL	36	19.8 mm	5.78 x 5.7 x 4.72 in (14.7 x 14.5 x 12.0 cm)	10/pk	Red	<input checked="" type="checkbox"/>
<b>HS120384</b>	15 mL	36	19.8 mm	5.78 x 5.7 x 4.72 in (14.7 x 14.5 x 12.0 cm)	10/pk	Purple	<input checked="" type="checkbox"/>
<b>HS120366</b>	50 mL	16	31.6 mm	5.78 x 5.7 x 4.72 in (14.7 x 14.5 x 12.0 cm)	10/pk	Natural	<input type="checkbox"/>
<b>HS120367</b>	50 mL	16	31.6 mm	5.78 x 5.7 x 4.72 in (14.7 x 14.5 x 12.0 cm)	10/pk	Blue	<input checked="" type="checkbox"/>
<b>HS120368</b>	50 mL	16	31.6 mm	5.78 x 5.7 x 4.72 in (14.7 x 14.5 x 12.0 cm)	10/pk	Red	<input checked="" type="checkbox"/>
<b>HS120369</b>	50 mL	16	31.6 mm	5.78 x 5.7 x 4.72 in (14.7 x 14.5 x 12.0 cm)	10/pk	Purple	<input checked="" type="checkbox"/>



Thin Polypropylene film 1.5-2 mL

Part No.	Tube	Wells	Well Dia.	L x W x H	Pkg.	Color	
<b>HS120370</b>	1.5/2.0 mL	81	12.1 mm	4.95 x 4.92 x 1.92 in (12.5 x 12.4 x 4.9 cm)	10/pk	White	<input type="checkbox"/>
<b>HS120371</b>	1.5/2.0 mL	81	12.1 mm	4.95 x 4.92 x 1.92 in (12.5 x 12.4 x 4.9 cm)	10/pk	Blue	<input checked="" type="checkbox"/>
<b>HS120372</b>	1.5/2.0 mL	81	12.1 mm	4.95 x 4.92 x 1.92 in (12.5 x 12.4 x 4.9 cm)	10/pk	Lt. Green	<input checked="" type="checkbox"/>
<b>HS120373</b>	1.5/2.0 mL	81	12.1 mm	4.95 x 4.92 x 1.92 in (12.5 x 12.4 x 4.9 cm)	10/pk	Purple	<input checked="" type="checkbox"/>



When finished with sample manipulation and ready for cold temperature storage, true north® freezer racks offer a variety of options to safely store samples in your freezer. You can choose from a wide selection of corrosion resistant stainless steel racks that are compatible with major brands or upright and chest freezers.

- All have handles for easy use
- Racks store boxes up to a 135 x 135 mm footprint and accommodate 2" tall boxes (up to 53 mm) as well as 3" tall boxes (up to 78 mm) depending on the rack mod

