

# MAGLab Series™

## Create lab storage space — off your bench top

With lab bench top space at a premium, common lab products such as pipettors, test tubes, Kimwipes®\* and miscellaneous accessories can be easily organized with Heathrow Scientific's magnetic lab storage system, the MagLab. MagLab includes four magnetic pieces: "MagRack" that holds (9) 15ml and (4) 50ml tubes; "Mag Wipe" that holds a standard size (11 cm x 21cm) Kimwipes; box; "MagPette" that holds two pipettors (universal design fits most common brands); and "MagBox" that holds miscellaneous items such as pens/pencils, notes, etc. Each rack is made of robust lightweight ABS plastic in bright and fun translucent colors and features four powerful magnets for a secure hold onto most metal surfaces or create your own holding point with the metal plate included for each segment. The MagLab Series is available individually or as a set of four.

Part No.	Description	L x W x H		Color
HS23500	maglab series, 4/pk	20.5 x 15.4 x 4 in	52.2 x 39.3 x 10.2 mm	■ ■ ■ ■ *
HS23501	magrack	5.1 x 4.4 x 2.5 in	13.1 x 11.2 x 6.4 mm	Teal ■
HS23502	magwipe	5.1 x 3.5 x 3.8 in	13.0 x 9.1 x 9.8 mm	Green ■
HS23503	magpette	5.1 x 4.4 x 2.9 in	13.1 x 11.2 x 7.6 mm	Purple ■
HS23504	magbox	5.1 x 3 x 4 in	13.0 x 7.8 x 10.2 mm	Orange ■

\*Assorted colors: Teal, Green, Purple and Orange

Kimwipes® is a registered trademark of Kimberly-Clark Corp.



Includes self-adhesive metal plates that adhere to any surface

Heathrow Scientific®

## Glove Box Holders

Keep gloves organized and ready for immediate use

- Free up bench top space by mounting holders on wall for easy dispensing of gloves
- Universal fit for most major brands of glove boxes
- Mounting screws included

Part No.	Description	Material	W x H x D
<b>HS23457</b>	1-box holder	Wire	140 x 209 x 108 mm
<b>HS23458</b>	3-box holder	Wire	250 x 454 x 108 mm

- Epoxy-coated steel wire racks are an economical solution



## Clearly Safe® Workstation Storage Bin

The essential combination accessory holder

- Organize your most commonly used laboratory disposables in one central location
- Acrylic unit has room to store gloves, Kimwipes®, and Parafilm®
- Use bins to store tubes, tips, etc.

Part No.	L x W x H
<b>HS234523</b>	265 x 137 x 440 mm

Parafilm® is a registered trademark of Bemis Inc.

## Pipette Workstation

Organize your pipette and tips

Keep pipette and disposable tips together in one convenient stand.

- Acrylic rack holds up to four pipettors on top shelf
- Dual bins house two different sizes of tips. Self-closing lids provide dust-free storage
- Oversized lip on lid allows for easy access even when wearing gloves

Part No.	L x W x H
<b>HS234631</b>	244 x 165 x 246 mm

