





**Store in ultra-low temperatures down to
-80°C with confidence**

True North® Slide Boxes offer a unique solution for storing slides down to -80°C. Durable polycarbonate construction with a rugged stainless steel clasp provides safe storage and transport for 3x1" (25x75mm) slides. The cryo slide box provides a 100 or 25-place capacity with uniform grooved slide slots to separate and prevent contact.

A permanent identification sheet is included in the hinged lid that corresponds to numbering on the temperature-resistant foam lining in the base. These features provide ideal inventory maintenance ability. The boxes are stackable.

Part No.	Capacity	Hinge	L x W x H	Color
HS15988H	100 slides	Nickle	208 x 175 x 34 mm	Purple 
HS120316	25 slides	Friction fit lid	141 x 88 x 35 mm	Blue 



Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (25 x 75 mm and 1 x 3 inch)

- Molded, rigid boxes securely hold slides in grooved slots
- Durable ABS plastic boxes are available with either cork or foam lining
- Both the lining and the slots are numbered to correspond to the slide inventory sheet



25 Place Slide Box

Easy access to slides with friction-fit cover

Part No.	Capacity	Lining	Hinge	L x W x H	Color
HS15989A	25 slides	Cork	Friction fit	141 x 88 x 35 mm	Blue <input checked="" type="checkbox"/>
HS15989B	25 slides	Cork	Friction fit	141 x 88 x 35 mm	Green <input checked="" type="checkbox"/>
HS15989C	25 slides	Cork	Friction fit	141 x 88 x 35 mm	Red <input checked="" type="checkbox"/>
HS15989E	25 slides	Cork	Friction fit	141 x 88 x 35 mm	White <input type="checkbox"/>
HS15989M	25 slides	Foam	Friction fit	141 x 88 x 35 mm	Blue <input checked="" type="checkbox"/>
HS15989P	25 slides	Foam	Friction fit	141 x 88 x 35 mm	Red <input checked="" type="checkbox"/>
HS15989R	25 slides	Foam	Friction fit	141 x 88 x 35 mm	White <input type="checkbox"/>



50 Place Slide Box

Thumb-latch lock mechanism

Available in cork lining only

Part No.	Capacity	Lining	Hinge	L x W x H	Color
HS15996A	50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	Blue <input checked="" type="checkbox"/>
HS15996B	50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	Red <input checked="" type="checkbox"/>
HS15996C	50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	White <input type="checkbox"/>



100 Place Slide Box

Closes securely with a rust-resistant, nickel-plated clasp and pin hinges

Part No.	Capacity	Lining	Hinge	L x W x H	Color
HS15994A	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Blue <input checked="" type="checkbox"/>
HS15994B	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Green <input checked="" type="checkbox"/>
HS15994C	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Red <input checked="" type="checkbox"/>
HS15994D	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Yellow <input checked="" type="checkbox"/>
HS15994E	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	White <input type="checkbox"/>
HS15994F	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Gray <input checked="" type="checkbox"/>
HS15994G	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Black <input checked="" type="checkbox"/>
HS15994M	100 slides	Foam	Nickle clasp	208 x 175 x 34 mm	Blue <input checked="" type="checkbox"/>
HS15994P	100 slides	Foam	Nickle clasp	208 x 175 x 34 mm	Red <input checked="" type="checkbox"/>
HS15994R	100 slides	Foam	Nickle clasp	208 x 175 x 34 mm	White <input type="checkbox"/>



100 Place Tall Slide Box

Stackable with 100-place premium slide boxes

Available in cork lining only

Additional space for larger 3 x 2 in

Part No.	Capacity	Lining	Hinge	L x W x H	Color
HS15979A	100 slides	Cork	Nickle clasp	208 x 175 x 60 mm	Blue <input checked="" type="checkbox"/>

Standard Microscope Slide Boxes

Constructed of durable polypropylene

- Slide boxes hold either 25 or 100 standard slides (25 x 75 mm and 1 x 3 inches)
- Hinged lids are held securely closed with a latch lock
- Stacking feet maximize stability
- Slide identification card
- Unlined and stackable



Part No.	Capacity	Hinge	L x W x H	Color
HS15990A	25 slides	Latch lock	141 x 92 x 30 mm	Blue <input checked="" type="checkbox"/>
HS15990B	25 slides	Latch lock	141 x 92 x 30 mm	Red <input checked="" type="checkbox"/>
HS15990C	25 slides	Latch lock	141 x 92 x 30 mm	White <input type="checkbox"/>
HS15991A	100 slides	Latch lock	210 x 169 x 37 mm	Blue <input checked="" type="checkbox"/>
HS15991B	100 slides	Latch lock	210 x 169 x 37 mm	Red <input checked="" type="checkbox"/>
HS15991C	100 slides	Latch lock	210 x 169 x 37 mm	White <input type="checkbox"/>