

## 4 Way Timer with Whiteboard and Pen

Keep track of time and important information with the 4 Channel Whiteboard Timer

### Features:

- Lays flat, stands up or use the magnetic back to attach to a metal surface
- Built in whiteboard area for writing
- Dry erase marker with eraser tip (included) helps keep track of crucial information
- Audible alarm for each channel
- Clock features count up or down to 24 hours 59 minutes
- Programs: hours/minutes; Displays: hours/minutes/seconds
- Up to 4 group timers can be set at the same time or used independently
- Last time set memory
- Key press confirmation feature gives an audible indication of data entry
- Sounds when time expires
- Batteries included (2 AAA)



Item No.	No. of Channels	Timing Capacity	Count Direction	Display Digits, H mm	Battery Included	L x W x H in	L x W x H cm
120365	Four	24 hr.	Up/Down	10	Two AAA	4.98 x 3.35 x .97	12.3 x 8.6 x 2.4

## Four Channel Timer/Clock

- Independently or simultaneously set time on four channels
- Audible and visual alarms signify elapsed time
- Easily toggle between channels and clock with top buttons
- Includes spring clip with magnetic back



Item No.	No. of Channels	Timing Capacity	Count Direction	Display Digits, H mm	Battery Included	L x W x H in	L x W x H cm
HS24670	Four	99 hr.	Down	14	Two LR44 1.5V	2.9 x 2.3 x 0.5	7.3 x 5.9 x 1.3

## Pocket Timer/ Stopwatch/Clock

- Four independent timing channels with visual and audible alarms
- Use as a stopwatch, a count up/down timer, or a clock
- Includes spring clip and magnetic back



Item No.	No. of Channels	Timing Capacity	Count Direction	Display Digits, H mm	Battery Included	L x W x H in	L x W x H cm
HS24490D	Four	100 hr.	Up/Down	6	One LR44 1.5V	2.4 x 2.8 x 1	6.1 x 7.1 x 2.5

## Triple Timer/Clock

- Three channels to count up/down with
- Simultaneous display, time-out, and 60 second audible alarm
- Channels may be started independently or simultaneously
- Timer automatically counts up after counting down to zero
- Memory recall



Item No.	No. of Channels	Timing Capacity	Count Direction	Display Digits, H mm	Battery Included	L x W x H in	L x W x H cm
HS24490M	Three	20 hr.	Up/Down	10.5	One AAA	3.4 x 2.6 x 0.9	8.6 x 6.6 x 2.2

## Multi-function Stopwatch/Clock

- Quick start, quick stop timer features time, split, hourly chime, calendar, alarm
- 100 cm neck cord



Item No.	No. of Channels	Timing Capacity	Count Direction	Display Digits, H mm	Battery Included	L x W x H in	L x W x H cm
HS24645	One	24 hr.	Up	8.5	One LR44 1.5V	3.4 x 3.1 x 0.8	8.7 x 7.8 x 2.1