

lab alert®

4 Way Timer with Whiteboard and Pen

Keep track of time and important information with the 4 Channel Whiteboard Timer.

Features:

- Lays flat, stands up or use the magnetic back to attach to a metal surface
- Built in whiteboard area for writing
- Dry erase marker with eraser tip (included) helps keep track of crucial information
- Audible alarm for each channel
- Clock features count up or down to 24 hours 59 minutes
- Programs: hours/minutes; Displays: hours/minutes/seconds
- Up to 4 group timers can be set at the same time or used independently
- Last time set memory
- Key press confirmation feature gives an audible indication of data entry
- Sounds when time expires
- Batteries included (2 AAA)

Part No.	No. of Channels	Timing Capacity	Count Direction	Display Digits, H	Battery	L x H x D
HS120365	Four	24 hr.	Up/Down	10 mm	Two AAA	4.98 x 3.35 x .97 in. 123 x 86 x 24 mm

Specifications

Resolution:	1/100 sec.
Clock:	yes
Alarm:	yes





a. Pocket Timer/Stopwatch/Clock

- Four independent timing channels with visual and audible alarms
- Use as a stopwatch, a count up/down timer, or a clock
- Includes spring clip and magnetic back



b. Triple Timer/Clock

- Three channels to count up/down
- Simultaneous display, time-out, and 60 second audible alarm
- Channels may be started independently or simultaneously
- Timer automatically counts up after counting down to zero
- Memory recall



c. Big-digit Dual Timer/Clock

- Display clock and timer simultaneously, or two independent timers
- Use one timer to count up, or use two timers to count down independently
- Audible and visual alarm



d. Multi-function Stopwatch/Clock

- Quick start, quick stop timer features time, split, hourly chime, calendar, alarm
- 100 cm neck cord



e. Four Channel Timer/Clock

- Independently or simultaneously set time on four channels
- Audible and visual alarms signify elapsed time
- Easily toggle between channels and clock with top buttons
- Includes spring clip with magnetic back

Part No.	No. of Channels	Timing Capacity	Count Direction	Display Digits, H	Battery	L x H x D
a HS24490D	Four	100 hr.	Up/Down	6 mm	One LR44 1.5V	61 x 71 x 25 mm
b HS24490M	Three	20 hr.	Down	10.5 mm	One AAA	86 x 66 x 22 mm
c HS24780	Two	24 hr.	Up/Down	25.5 mm	One AAA	142 x 97 x 20 mm
d HS24645	One	60 min.	Up	8.5 mm	One LR44 1.5V	87 x 78 x 21 mm
e HS24670	Four	99 hr.	Down	14 mm	Two LR44	73 x 59 x 13 mm

Specifications	
Resolution:	1/100 sec.
Clock:	yes
Alarm:	yes