

## Wash Bottles

### Color-coded for easy identification

The closure and stem are molded in one piece making them leak-proof

The angled neck and tip provide a uniform stream of liquid and the tip can be cut back to increase the flow

- Wide bottle mouth makes filling easier
- Bottle is made from LDPE with a PP closure

Pack of 6

Item No.	Cap/Lid Color	Capacity
120245	White <input type="checkbox"/>	500 mL
120246	Red <input type="checkbox"/>	
120247	Yellow <input type="checkbox"/>	
120248	Green <input type="checkbox"/>	
120249	Blue <input type="checkbox"/>	



## Labeled Wash Bottles

### Clear and informative product safety information

Wash bottles safety labeled and self-venting

The information provided on the front of each bottle is specific to the product contained and includes the following — the Chemical Abstract Service (CAS) No., National Fire Protection Association (NFPA) codes and symbols, Protective Equipment Symbol (PES), Target Organ Information and Hazard Codes

There is also the option of using other liquids within the bottles and creating the appropriate safety information on the outside of the bottle by using the 'self-labeling' versions

- Wide bottle mouth makes filling easier
- Sodium Hypochlorite bottle made from natural LDPE for light protection and has a HDPE closure
- All other bottles are made from LDPE with a PP closure
- Non-metallic vent that prevents the risk of solvent leakage
- Six named chemicals, a self-labeling\*\* and assorted pack of bottles to choose from

Item No.	Cap/Lid Color	Description	Capacity
120250	Red	Acetone	500 mL
120251	Clear	Distilled Water	
120252	Yellow	Isopropanol	
120253	White	Ethanol	
120254	Green	Methanol	
120255	Blue	Sodium Hypochlorite	
120256	White	Self-labeling**	
120257	Assorted*	Assorted*	

\* Assorted Pack includes one of each: acetone, distilled water, isopropanol, ethanol, sodium hypochlorite & a self-labeling bottle

\*\*Self-labeling bottles must only be used with chemicals that are safe and compatible with LDPE. Labeling stickers included

